




Angela Hinchie, Elena Vialetto

Vivlion's PRCISR™ CRISPR platform


PRCISR™ CRISPR is Vivlion's **full-service** CRISPR-enabled discovery platform that solves key challenges of CRISPR screening, thereby optimizing the target identification process.



High Throughput CRISPR Screening Service

Since the inception of CRISPR-Cas editing, scientists have recognized the value in using it for forward genetic screens. Conducting a full-sized screen is a monumental task that can yield exciting genetic findings – but only if done correctly. Vivlion offers support for each step of the screening process, from initial design to the data analysis.

Diverse and uniformly distributed CRISPR libraries




What is uniformity and why does it matter?

Library skew measures the ratio of the most to the least represented gRNAs in a library. Unequal distribution and guide loss can hinder bioinformatic analysis and lead to lower confidence results. Vivlion's proprietary 3Cs synthesis method results in more uniform libraries with fewer dropped guides than conventional production methods. Conducting the same screen with a more uniform library yields more high confidence hits.

Alexandria: Choosing the best guides for the job


Vivlion's Alexandria library was created by analyzing empirical data on gRNA activity and identifying the most active guides over a wide range of contexts. This improved library is available in single and dual-targeting formats and further improves the CRISPR screening process.



Vivlion proudly works with: **REPAIR** THERAPEUTICS


Low Skews facilitate downscaling and parallelization

Conventionally cloned CRISPR libraries have high skews that lead to low confidence in potential hits as shown by Vivlion's partner Repair Therapeutics. The conventional library had a higher skew and showed fewer hits when screened at 100x as compared to 400x. The TKOv3-3Cs maintained a high hit retention rate even at a 100x as compared to the conventional library at 400x.



Dual-targeting increases editing efficiency up to 10x

Targeting the same gene with two gRNAs increases the probability of efficient gene knock out by a factor of 10x as shown by Vivlion's partner Repair Therapeutics.



Figures courtesy of Repair Therapeutics.

PRCISR™ CRISPR Services

Screens

- Experimental planning and design
- Combinatorial options
- Multiple readouts supported
- High throughput lab capacity
- Experienced CRISPR scientists

Libraries

- PCR-free production
- Highly uniform
- Customizable
- Validated gRNAs
- Available in single, dual and combinatorial formats

Bioinformatics

- *In-silico* library design
- Customizable for singular and combinatorial edits
- Established pipelines
- Genetic interaction maps
- Data integration

References: Wegner et al., eLife 2019 | Henkel et al., BMC Biology, 2020 | Diehl et al., NAR, 2021 | Wegner et al., Bioinformatics, 2023